



製品仕様

●材質

- ・バレル : 環状ポリオレフィン (COP)
- ・チップキャップ : ブチルゴム
- ・ピストン : ブチルゴム、フッ素ラミネート
- ・プランジャーロッド : ポリカーボネート、ポリプロピレン、ポリスチレン
- ・グリップ : ポリプロピレン

●バレルサイズ表

<ルアーロックタイプ>

| 容量 (mL) | 外径 (mm) | 内径 (mm) | 全長 (mm) |
|---------|---------|---------|---------|
| 0.5 | 9.40 | 4.65 | 64.4 |
| 1.0 | 9.40 | 6.35 | 64.4 |
| | 11.20 | 6.35 | 64.4 |
| 2.25 | 11.20 | 8.65 | 64.4 |
| 3.0 | 11.20 | 8.65 | 84.4 |
| 5.0 | 15.05 | 12.45 | 79.0 |

<ルアーコーンタイプ>

| 容量 (mL) | 外径 (mm) | 内径 (mm) | 全長 (mm) |
|---------|---------|---------|---------|
| 1.0 | 9.15 | 6.35 | 64.5 |
| 2.25 | 10.85 | 8.65 | 64.4 |
| 5.0 | 15.05 | 12.45 | 79.0 |

ClearJect® WN

【認証番号】 225AGBZX00025000


容量 : 1mL、2.25mL

針 : 27G、29G (½" Thin wall)
22G-26Gも対応可



製造工程

クリーンな環境での全自動ライン製造



包装仕様

- TASPAC® (RTUタイプ)
- γ線滅菌済み



※ご要望に応じた専用仕様にも対応致します。



Product Specification

● Materials

- Barrel : Cyclo Olefin Polymer
- Tip Cap : Butyl Rubber
- Piston : Chlorinated Butyl Rubber
or Fluoro-Laminated Rubber
- Plunger Rod : Polypropylene, Polystyrene or Polycarbonate
- Grip : Polypropylene

● Barrel Size

<Luer Lock Type (LL)>

| Volume (mL) | O.D. (mm) | I.D. (mm) | Total Length (mm) |
|-------------|-----------|-----------|-------------------|
| 0.5 | 9.40 | 4.65 | 64.4 |
| 1.0 | 9.40 | 6.35 | 64.4 |
| | 11.20 | 6.35 | 64.4 |
| 2.25 | 11.20 | 8.65 | 64.4 |
| 3.0 | 11.20 | 8.65 | 84.4 |
| 5.0 | 15.05 | 12.45 | 79.0 |

<Luer Cone Type (LC) >

| Volume (mL) | O.D. (mm) | I.D. (mm) | Total Length (mm) |
|-------------|-----------|-----------|-------------------|
| 1.0 | 9.15 | 6.35 | 64.5 |
| 2.25 | 10.85 | 8.65 | 64.4 |
| 5.0 | 15.05 | 12.45 | 79.0 |

ClearJect® WN

[Approval No.225AGBZX00025000]


Volume : 1mL, 2.25mL

Needle : 27G, 29G (½" Thin wall)
also applicable with
22G, 23G, 25G and 26G



Production Process

Fully automated manufacturing under the clean condition.



Packing

- TASPAC®(RTU packing)
- Gamma Sterilization Finished



※Your required packing applicable.

If you are an overseas customer (outside Asia), please contact Gerresheimer as shown below:

URL : <https://www.gerresheimer.com/>

Phone : +49 5223 164 309 E-mail : clearject@gerresheimer.com

Person to contact : Mr. Bernd Zeiss