


プラスチック製 針付プレフィル用シリンジ PAT.

ClearJect® WN

【認証番号】225AGBZX00025000



	容量 (mL)	外径 (mm)	内径 (mm)	シリンジ長 (mm)	RNS 装着全長 (mm)	フランジ外径 (mm)
ClearJect® WN 1mL long	1.0	8.15	6.35	76.7	81.7	13.8
ClearJect® WN 2.25mL	2.25	10.85	8.65	79.3	84.3	17.75
ClearJect® WN 2.25mL SRF	2.25	10.85	8.65	79.3	84.3	14.7


包装仕様

- TASPAC® (RTU タイプ)
ガンマ線滅菌済み



特長

- 接着剤不使用
- 単回使用、使い捨て
- 割れ難い
- 高い水蒸気バリア
- 専用シリコーン採用による低い穿刺抵抗



Plastic Prefillable Syringe with Needle PAT.

ClearJect® WN

【 Approval No.225AGBZX00025000 】

Needle:

27G 1/2inch (Thin wall)

also applicable with

22G, 23G, 25G, 26G and 29G

Stainless (SUS304)

Insert molding (Glue-free)

Barrel:

Dimensions are in accordance
with ISO 11040-4, (ISO 11040-6)
or custom-made size.
Cyclo Olefine Polymer

RNS (Rigid Needle Shield):

Outer Cap : PP

Inner Cap : Chlorinated Butyl Rubber

2.25mL

1mL

	Volume(mL)	O.D.(mm)	I.D.(mm)	Syringe Length(mm)	Total Length with RNS(mm)	Flange O.D.(mm)
ClearJect® WN 1mL long	1.0	8.15	6.35	76.7	81.7	13.8
ClearJect® WN 2.25mL	2.25	10.85	8.65	79.3	84.3	17.75
ClearJect® WN 2.25mL SRF	2.25	10.85	8.65	79.3	84.3	14.7

Packaging

- TASPAC® (RTU Type)
- Gamma Sterilization Finished



Advantages

- Glue free
- Usable only once
- Hard-to-break
- High barrier property of water vapour
- Lower value of Needle Puncture Resistance by exclusive use silicone

