

VF-SIN Syringe

ガラス製針付プレフィル用シリンジ



製品仕様

製品名	容量 (mL)	全高 (mm)	外径 (mm)	内径 (mm)	針長 (in)
VF-SIN Syringe 0.5mL	0.5	47.6	6.85	4.65	1/2
VF-SIN Syringe 1mL std.	1.0	35.7	10.85	8.65	1/2 , 5/8
VF-SIN Syringe 1mL long	1.0	54.0	8.15	6.35	1/2
VF-SIN Syringe 1.5mL	1.5	43.2	10.85	8.65	1/2 , 5/8
VF-SIN Syringe 2.25mL	2.25	54.4	10.85	8.65	1/2

● 材質

パレル: ガラス

針: ステンレス(Thin wall)

ガス透過キャップ

: イソプレンゴム

: 熱可塑性エラストマ(TERNS)

包装仕様

● RTF[®] (Ready-to-Fill)

滅菌済み

RTF[®]は Gerresheimer Buende 社の登録商標です。



VF-SIN Syringe

Glass Prefillable Syringe with Needle



Product Finish

Product name	Capacity (mL)	Height (mm)	Outer Diameter (mm)	Inner Diameter (mm)	Needle Length (in)
VF-SIN Syringe 0.5mL	0.5	47.6	6.85	4.65	1/2 , 5/8
VF-SIN Syringe 1mL std.	1.0	35.7	10.85	8.65	5/8 , 1
VF-SIN Syringe 1mL long	1.0	54.0	8.15	6.35	1/2 , 5/8
VF-SIN Syringe 1.5mL	1.5	43.2	10.85	8.65	1/2 , 5/8
VF-SIN Syringe 2.25mL	2.25	54.4	10.85	8.65	1/2

- Material
- Barrel: Glass
- Needle: Stainless Steel (Thin wall)
- Gas Permeable Cap : Isoprene Rubber
- : Thermo Plastic Elastomer (TERNS)

Packaging specification

●RTF[®] (Ready-to-Fill)

Sterilized

RTF[®] is a registered trademark of Gerresheimer Buende GmbH.



If you are an overseas customer (outside Japan), please contact Gerresheimer as shown under:

URL: <https://www.gerresheimer.com/>

(Europe & Asia) Phone: +49 5223 164 239

(Americas) Phone: +1 856 405 9737

E-mail: info-syringes@gerresheimer.com

E-mail: info-syringes@gerresheimer.com

Person to contact : Ms. Beate Desjardins

Person to contact : Ms. Karen Kanyok

